



UPDATED
VERSION
2026

CASIA2

USER EXPERIENCE

You + eye. We care.

"SEEING CASIA2 USED IN
REAL CLINICAL SETTINGS
IS AMAZING - ITS WIDE
RANGE OF APPLICATIONS
AND EASE OF USE ARE
TRULY IMPRESSIVE."

Udo Schmidt

HEAD OF PRODUCT MANAGEMENT/
R&D





+ *Table of contents*

CONTENT FROM CASIA2 USERS

Felix Gonzalez-Lopez MD, PhD	p. 6	Dr. Teresa Sandmayr, FEBO /	p. 42
Dr. Farzona Abdurakhimova /	p. 12	Dr. Aron Roman Czappek	
Dr. Khabibullaeva Nigora Khusanovna		Dr. Stephan Szegedi	p. 44
Prof. Dr. Timo Eppig / Dr. Eckhard Becker /	p. 16	Ana Trpeska Boshkoska MD	p. 48
Violetta Müller		Dr. Jorge Vila Arteaga	p. 50
Dr. Luca Gualdi	p. 18	Prof. Adam Wylęgała	p. 52
Dr. Carlos Gutiérrez	p. 20	Azamat Farkhadovich Yusupov, MD, DSc, Prof /	p. 56
Dr. Merieme Harouch	p. 22	Ikbol Fazilovich Saliev, MD, PhD /	
Priv. Doz. Dr. Nino Hirnschall, PhD, FEBO, MhBA /	p. 24	Farzona Abdurakhimova, MD	
Christopher Reischl BSc.			
Dr. Clemens Kaindlstorfer, FEBO /	p. 26		
Dr. Aron Roman Czappek			
William Kay / Ayan Abdullahi / Madyan Qureshi	p. 28		
Dr. Sophia Reifeltshammer, MD /	p. 30		
Dr. Klaus Straßmair, MD			
Ignacio Rodriguez-Una MD, PhD	p. 34		
Dr. Pau Romera	p. 38		

CONTENT FROM TOMEY

CASIA2 features	p. 4
Glaucoma application	p. 37
ICL surgery	p. 41
Join this project	p. 47
Artificial prosthetic eyes	p. 59

CASIA2

Cornea/Anterior Segment OCT

Optimise clinical decision-making with the CASIA2 technology

- + Testing application for cataract/glaucoma/cornea surgery
- + Special glaucoma analysis (STAR 360°)
- + Advanced imaging with deep scanning depth (13 mm)
- + Very fast scanning speed (50,000 A-scans/second)
- + Corneal topography, trend analysis, IOL choice, and calculation
- + Lens shape analysis incl. tilt and decentration
- + Phakic IOL simulation
- + ICL Pre- and Post-OP application
- + Colour fixation camera

Special thanks
to all the experts
who submitted their
practical cases.

In this brochure, we proudly present a collection of diverse patient experiences that showcase the versatility and clinical utility of CASIA2 Cornea / Anterior Segment OCT. Dive into the insights, outcomes, and lessons learned from these cases, and see firsthand how CASIA2 is transforming ophthalmic care.





Felix Gonzalez-Lopez MD, PhD


Medical Director,
Clinica MIRANZA, Madrid (Spain)

Introduction

1. Anterior Segment Optical Coherence Tomography (AS-OCT) represents an unparalleled tool for assessing the anterior pole of the eye, providing a novel approach to studying the anterior segment. This innovation has led to a significant paradigm shift in preoperative calculations for sizing posterior chamber phakic lenses.


One of the latest enhancements to the CASIA2 is the **Post-op ICL** application that facilitates quick and straightforward postoperative assessment of posterior chamber phakic lenses. I have had the privilege of actively contributing to its development. With this application, we can now not only visualise but also measure the dynamic vault of the ICL, observing how it varies with external lighting conditions. Moreover, it enables graphical assessment of the relationship between pupil size, the optic zone of the phakic lens, and the status of the iridocorneal angle. The imagery it provides under scotopic and photopic conditions, easily exportable to medical records, fulfills the criteria every clinician seeks when examining an eye implanted with an ICL.

With the introduction of this new application, the **CASIA2** unquestionably emerges as the most potent imaging device available for postoperative evaluation of ICLs.




Low central vault

2. In 2018, we introduced what we termed the **dynamic concepts of vault** (Gonzalez-Lopez F, et al. Dynamic assessment of light-induced vaulting changes of implantable collamer lens with central port by swept-source OCT: Pilot study. Trans Vis Sci Tech. 2018;7(3):4. <https://doi.org/10.1167/tvst.7.3.4>). In this paper, we outlined the novel parameters **vault interval (VI)**, which we defined based on central vault values during maximum mydriasis and maximum miosis following light-induced changes in pupil diameter, and **vault range (VR)**, calculated as the absolute difference between the VI values. In our initial pilot study, we found that the mean VR was $167 \pm 70 \mu\text{m}$. Through a collaborative effort with TOMEY, we were able to integrate the dynamic vault concept into an application, the **Post-op ICL**, now accessible in CASIA2.



3. The dynamic images captured by the **video application**, recorded under controlled scotopic and photopic ambient lighting conditions, offer us an authentic perspective on the interactions between the ICL and different structures of the anterior ocular segment. Combined with the **Post-op ICL** application's ability to quantify these changes, this leads us to introduce another new concept: the **relativism of the vault**. Consequently, it enables us to determine, on a case-by-case basis, whether pupillary dynamics, iridocorneal angles, and vault are optimal, or conversely, whether explantation and replacement with an ICL of a different size is necessary.






High vault






Low peripheral vault



Alternatively, as demonstrated in the example presented, rotating the phakic lens vertically can reduce vault and thereby improve pupillary dynamics. This surgical maneuver successfully alleviated the dysphotopsias experienced by the patient in this case.



4. Preoperative evaluation of **iris morphology**, its displacement under varying external lighting conditions, and **crystalline lens rise (CLR)** undoubtedly aids in selecting the appropriate ICL size. In this instance, a notably low CLR, coupled with a deep anterior chamber depth (ACD) and concave iris shape, will impact the choice of ICL size. It is noteworthy to observe the enhancement in pupillary dynamics and the modification in the relationship between the iris and the crystalline lens following ICL implantation.



Case:

CASIA2 in Cataract Surgery: From Preoperative Planning to Postoperative Confirmation

This case illustrates how CASIA2 can contribute to a high-quality cataract surgery workflow, not only by supporting preoperative assessment and IOL planning, but also by providing objective postoperative confirmation of the anatomical and refractive outcome.

The patient was a 60-year-old male with high myopia in both eyes and visually significant cataracts. Preoperative refraction in the right eye was $-8.00 / -2.00 \times 15^\circ$, with a corrected decimal Snellen visual acuity of 0.32. Slit-lamp examination showed a corticonuclear cataract with a dense posterior subcapsular component. Axial length was 28.18 mm.

In this type of eye, careful preoperative evaluation is particularly important. CASIA2 provided a detailed analysis of the cornea, anterior chamber, crystalline lens, and biometric parameters, helping to confirm the anatomical conditions for surgery and to refine the preoperative planning.



Felix Gonzalez-Lopez MD, PhD
Medical Director,
Clinica MIRANZA, Madrid (Spain)




Figure 1: Preoperative Total Analysis

Comprehensive preoperative CASIA2 Total Analysis of the anterior segment. The examination documents corneal morphology, pachymetry, anterior chamber parameters, pupil size, crystalline lens position, lens tilt, and decentration before cataract surgery.




Figure 2: Preoperative Cataract Planning and IOL Calculation

CASIA2 preoperative cataract planning screen showing axial length (data imported from the biometer), anterior chamber depth, lens thickness, white-to-white distance, keratometric data, and IOL calculation. These data supported surgical planning.

The patient underwent femtosecond laser-assisted cataract surgery (FLACS). A Tecnis PureSee (Johnson & Johnson) ZEN00V +9.00 D intraocular lens was implanted.

Postoperatively, the patient achieved an excellent functional and refractive result, with uncorrected visual acuity of 1.00 and plano refraction. The postoperative analysis confirmed a stable and well-positioned IOL, with minimal tilt and decentration, an expected increase in anterior chamber depth after crystalline lens removal, and stable corneal parameters.



Figure 3: Postoperative Cataract Analysis

Postoperative CASIA2 cataract analysis comparing the preoperative and postoperative status. The examination confirms the anatomical changes after cataract surgery, including deepening of the anterior chamber and preservation of corneal stability.



Figure 4: Postoperative IOL Analysis

CASIA2 postoperative IOL analysis showing the position of the implanted intraocular lens. The IOL was well positioned, with minimal tilt and decentration, consistent with the excellent refractive and visual outcome. Further information regarding CASIA2-based quantification of crystalline lens and IOL tilt and decentration can be found in: González-López F, Delgado-Tirado S, Pérez-Izquierdo R, González-González D, Medel D, Sáez-Corrales A, Péláez-Gutiérrez EC. Tilt and Decentration of the Crystalline Lens and Intraocular Lens After Femtosecond Laser-Assisted Cataract Surgery. *J Cataract Refract Surg.* 2026;52(2):140-147. doi:10.1097/j.jcrs.0000000000001770.

In this case, CASIA2 was not required to solve a complex complication, but rather to support excellence in a routine cataract case. It provided a complete, objective, and visually intuitive documentation of the eye before and after surgery. This type of information is valuable for preoperative planning, postoperative quality control, patient documentation, and clinical confidence.

The final outcome was fully consistent with the surgical plan: excellent uncorrected visual acuity, plano refraction, stable corneal morphology, and a well-centred intraocular lens.

Case:

Post-ACXL Corneal Edema in Advanced Keratoconus


Patient: 27-year-old male
 Diagnosis: Bilateral keratoconus — OD: stage III, OS: stage IV
 Procedure: Accelerated corneal cross-linking (ACXL),
 Epi-Off (9 mW/cm², 10 minutes) on OD

Clinical Context


Following ACXL, the patient presented with visible stromal edema and Descemet's membrane folds. Slit-lamp visualisation was limited by corneal haze. The need for accurate, reproducible, and non-invasive imaging was clear. CASIA2 provided the precision and depth of information necessary to track the healing process over time.



Farzona Abdurakhimova, MD
 Ophthalmic surgeon and researcher,
 SAIF Optima Eye Clinic,
 Tashkent (Uzbekistan)




Dr. Khabibullaeva Nigora Khusanovna
 Ophthalmic surgeon and researcher,
 SAIF Optima Eye Clinic,
 Tashkent (Uzbekistan)




CASIA2 scores before Epi-Off corneal ACXL

Monitoring with CASIA2

	CASIA2 Findings	Clinical Relevance
Day 1	Stromal edema, Descemet's folds	Initial post-op response
Day 7	Edema unchanged	Safe observation confirmed



CASIA2 scores on day 1 after Epi-Off corneal ACXL



CASIA2 scores on day 7


Why CASIA2?

The CASIA2 swept-source OCT is a transformative advance in anterior segment imaging. With its unparalleled resolution, deep scanning depth, and ultra-fast acquisition, it offers a level of clarity, confidence, and control that is simply unmatched.


Whether it's for keratoconus diagnostics, angle assessment, ICL planning, or post-op monitoring, CASIA2 gives clinicians the tools they need to see more, understand in more depth, and act sooner.

Monitoring with CASIA2


	CASIA2 Findings	Clinical Relevance
Day 14	Folds resolved, edema decreasing	Signs of spontaneous recovery
Day 30	Clear cornea, regular epithelium	Clinical resolution
2 Months	Stable morphology, VA 0.3	No complications
9 Months	VA 0.5-0.8, normal topography	Long-term stabilisation





CASIA2 scores on day 14



CASIA2 scores on day 30



CASIA2 scores 2 months after Epi-Off corneal ACXL



CASIA2 scores 9 months after Epi-Off corneal ACXL



Key Contributions of CASIA2 in this Case

- + Enabled detailed layer-by-layer imaging, even in the presence of haze
- + Provided objective validation of natural recovery
- + Helped avoid unnecessary treatment through precise monitoring
- + Ensured safe long-term follow-up with clarity and consistency

Conclusion: More Than a Device – A New Standard of Care

The CASIA2 is not just an imaging device – it is a revolution in anterior segment diagnostics. In this case, it offered what no slit lamp or ultrasound could: a non-invasive, highly detailed, and dynamically trackable view of the cornea's healing journey.

Its ability to visualise microstructural changes with such depth and accuracy makes it indispensable in modern ophthalmology. For corneal surgeons, refractive specialists, and glaucoma experts alike, CASIA2 raises the standard, shortens decision-making time, and enhances patient outcomes.

Once you start working with CASIA2, you will no longer imagine practicing without it. It is more than a tool – it is a clinical ally, a precision instrument, and for many of us it is the eye behind our eyes.



Prof. Dr. Timo Eppig

Adjunct professor of Experimental Ophthalmology at Saarland University, Homburg (Germany), Senior Consultant and Head of Research and Development at AMIPLANT GmbH, Schnaittach (Germany)

Dr. Eckhard Becker

Senior ophthalmic surgeon and co-owner at Augentagesklinik Oranienburg MVZ GmbH, Oranienburg (Germany) – a private eye clinic with several subsidiaries nearby Berlin (Germany)

Violetta Müller

Resident at Augentagesklinik Oranienburg MVZ GmbH, Oranienburg (Germany)

Case:

CASIA2 for Measuring the Ciliary Sulcus in Order to Determine the Size of an Accommodating Intraocular Lens

Introduction

The Lumina accommodating intraocular lens is a dual-optic lens system intended for implantation in the ciliary sulcus. The mechanism of action is based on the Alvarez lens principle. The lens is directly driven by the ciliary muscle, which laterally displaces the two lens components. The performance of the Lumina is directly depending on a perfect fit in the ciliary sulcus and they are therefore custom-made for each eye. Therefore, precise measurement of the sulcus dimensions is crucial for the procedure's success.

Ultrasound imaging provides poor resolution, even with high-frequency probes, so instead we use the CASIA2 for precise direct measurement of the ciliary sulcus. This is not only more comfortable for the patient but also delivers superior image quality.

Methods

We created a set of Scan Types using up to 64 averaging frames to reduce image noise. This improved the image quality in the ciliary sulcus so that the measurement tool in the 2D Analysis mode can be used to directly measure the width of the ciliary sulcus (Figure 1). It was measured in three different meridians.




Figure 1: Preoperative measurement of the width of the ciliary sulcus and the distance of the sulcus plane to the corneal vertex

The Lumina is then manufactured using these measurements, and the accurate fit of the lens was assessed postoperatively using the same method; it showed a perfectly placed lens in the ciliary sulcus (Figure 2). The performance of the Lumina was then evaluated by the PowerRef3 (PlusOptix GmbH, Nürnberg). It exhibited an accommodative amplitude of up to 2.5 D with a stimulus of 3 D.





Figure 2: Postoperative image of the Lumina in the eye

Conclusion

The CASIA2 is the ideal instrument for pre- and postoperative assessment of the ciliary sulcus: It provides superior image quality compared to ultrasound imaging and is also faster and less burdensome to the patients.





Dr. Luca Gualdi

Ophthalmologist, Anterior Segment
Surgeon at Diagnostica Oftalmologica
e Microchirurgia Oculare, Rome (Italy)

Case 1:

Vogt Koianagi Harada


Anterior segment manifestation of a Vogt Koianagi Harada uveitis. The early diagnosis was made thanks to the AS-OCT CASIA2 image: a choroidal effusion moved towards the ciliary muscle creating a secondary acute glaucoma in a preexisting narrow angle. The patient started immediately the right therapy (high dose intravenous steroids for the panuveitis and acetazolamide tablets with a combination of brimonidine and timolol) and the case was solved in few days.



Case 2:

LAPR (Laser Presbyopia Reversal)

Erbium laser effects of radial cuts into the scleral tissue to reduce ocular bio-mechanics and rigidity to restore the accommodation in a presbyopic patient. The deepness of the scleral tissue is still visible almost 30 years after LAPR surgery. Thanks to the swept-source AS-OCT CASIA2 images, we can see that there was not a complete perforation until the blue choroid. Given that LAPR is still performed manually, complications may occur and the scars had a different deepness being done manually. This old and amazing experimental technique will soon be replaced by a new OCT and computer-guided, safer and more efficient called LSM (Laser Scleral Microporation).



From my practice – what makes CASIA2 stand out

I enjoy the CASIA2 because being a swept-source OCT you can go deeper and see better the whole anterior segment even until the posterior capsule of the lens.

This makes an easier way to stage all kinds of cataracts. The integration with the OA-2000 biometer can give an instant calculation with a lot of formulas, including also the pre- (and post-) tilting of the lens and the angle kappa which can be an inclusion or exclusion for a premium IOL. The software is always updated and you can play with the internal toric IOL calculator which includes also the posterior and total cylinder. You can also choose the right astigmatism axis from the fourier and the okulix formula which I found very effective, especially in irregular corneas (e.g. keratoconus, irregular astigmatism or post-refractive, including post-RK).

The complete corneal trend analysis also makes follow-up easier, especially in my keratoconus patients where also the epithelial map is included.

With the newest software you can also follow up the narrowing of an irido-corneal angle, which can suggest or not a prevention for an acute glaucoma with a laser iridotomy or a lens replacement.

My favourite software is the inclusion of the Nakamura and KS formulas which simplify the selection for a correct vaulting in posterior chamber phakic IOLs.

The CASIA2 acquisition is very quick, easy and repeatable and you can play also in 3D images and with movies which can be very useful, especially to detect the dynamic vaulting in cases where a posterior phakic IOL is with a low vaulting.



Dr. Carlos Gutiérrez

Medical Director,
Clínica Vista Gutiérrez Amorós,
A Coruña (Spain)

The use of anterior chamber OCT has increased its potential over the years, providing more and higher-quality information. The CASIA2 has allowed us to offer a greater guarantee in the results of glaucoma surgeries by observing the effectiveness of these procedures instantly.

Case 1:

Assessing ESNOPER V2000 Implant Permeability

Glaucoma implants serve to drain the aqueous humor in order to reduce intraocular pressure. For this reason, it is important to verify that these implants are permeable and effective.



As seen in the images, CASIA2 system allows us to precisely evaluate the permeability of the ESNOPER V2000 scleral implant. Using the “bleb” measurement system the valve’s permeability in all the sectors can be evaluated, and the intrascleral space can be analysed to determine if it is large enough.



Case 2:

Evaluating Preserflo® Implant Efficiency

In devices used for the treatment of glaucoma, such as the Preserflo® implant, it is important to verify its position, depth, and effectiveness in reducing intraocular pressure through filtration. The use of the CASIA2 allows us to evaluate these parameters with a single measurement and a single image. With this, we can verify from the end of the implant in the anterior chamber to its functioning at the scleral level and the filtration capacity that is being produced. This undoubtedly provides a greater guarantee in the results of the surgery in the immediate postoperative period and in the control of the pathology in the long term.



Case:

The First Asymmetric Ring AS6-300/150° Implantation in Keratoconus

Refraction topo (axial power map) BCVA: 5/10 with (-3.50; 65)


Refraction on TOPO (Fourier Index asymmetry)

Refrac: BCVA (-3.50 [A]85°) = 1 elt of 8/10



Pr. As. Merieme Harouch

Cornea and refractive surgeon,
Casablanca cornée kératocône CCK,
Casablanca Laser Vision CLV,
Casablanca (Morocco)



Asymmetric corneal ring segment

Apical diameter: 6,80 mm

Base width: 800µ

Arc length: 150°

Thickness: 250-350-250

Tunnel parameters

Interior diameter: 5,95 mm


Exterior diameter: 7,35 mm



Interior depth = exterior depth:

80% thickness at implantation site

Comparative table of VA and Kmax

	UCVA	BCVA	Kmax
<i>Before ICR</i>	4/10	(-3.50 [A]85°) 1 elt of 7/10	51.30
<i>After ICR</i>	3 elts of 7/10	(-1.00 [A]115°) 2 elts of 8/10	48.93





My conclusions

- + The Asymmetric ring AS6-150° 250/350/250 is effective in flattening the cornea.
- + The flattening is more important at the top of the ectasia where the ring is thickest (350).
- + The flattening gradually decreases at the periphery of the ectasia, where the two ends of the ring are thin (250μ).
- + The centering of the asymmetric ring on the axis of the Fourier axial asymmetry gave a good functional and topographic results.

Effectiveness on axial asymmetry and therefore on aberrations (Coma), in addition to effectiveness on the SE.

Case 1:

Cyst at the Limbus




**Priv. Doz. Dr. Nino Hirschall, PhD
FEBO, MhBA**

Department of Ophthalmology,
Kepler University Clinic, Linz (Austria)





Christopher Reischl BSc.

Department of Ophthalmology,
Kepler University Clinic, Linz (Austria)




In this case of stem cell deficiency, cysts occurred on both eyes at the area of the limbus.



With the CASIA2 it was possible to observe corneal changes and measure the size and shape of the cysts from different angles. This resulted in a clearly visualised picture of the size of the cysts and enhanced the depth of information. Therefore, it was possible to evaluate the size and the lack of infiltration into the surrounding tissue, and then get a better understanding of the content of the cyst.

Case 2:

Repeatedly Operated Eye





Dr. Clemens Kaindlstorfer, FEBO

Department of Ophthalmology,
Barmherzige Brüder,
Konventhospital, Linz (Austria)





Dr. Aron Roman Czapek

Department of Ophthalmology,
Barmherzige Brüder,
Konventhospital, Linz (Austria)

Case:



Postoperative Angle Closure After Pars-Plana Vitrectomy in Pseudophakic Eye

This case involves a 71-year-old male patient who underwent pars-plana vitrectomy due to a retinal detachment in his left eye. The patient was already pseudophakic at the time of the surgery. Additionally, a YAG capsulotomy had already been performed. After an initially uneventful postoperative course, a complete angle closure developed on the second postoperative day due to a gas bubble in the posterior chamber, between the anterior IOL surface and the posterior iris, as a result of an absent capsular barrier. The intraocular pressure was 50 mmHg.



Figures 1 and 2: CASIA2 images of the left eye with the gas bubble in the posterior chamber leading to a 360° pupillary block

Subsequently, maximal local pressure-lowering therapy was administered, along with intravenous mannitol and acetazolamide. Furthermore the patient was instructed to adopt a prone position. Additionally, pilocarpine was applied to constrict the pupil.



Figures 3 and 4: Images of the same patient after conservative treatment and shifting of the gas bubble into the anterior chamber

The following day, the angle closure had resolved, and the gas bubble had shifted to the anterior chamber. The intraocular pressure was within the normal range at 11 mmHg. This case illustrates the possibility of avoiding a potentially difficult surgery through sophisticated conservative methods.

Case 1:

Descemet's Membrane Detachment (Pre-Op and Post-Op)

Patient had complicated cataract surgery with a small Descemet's membrane scroll intraoperatively. The patient was seen two days postoperatively and the cornea was very hazy, making it difficult to assess. The CASIA2 OCT demonstrated a clear Descemet's membrane detachment. The patient was brought back to theatre the same day and had an air tamponade injected to re-attach the Descemet's membrane. They were followed up with subsequent CASIA2 OCT scans to demonstrate successful re-attachment.



William Kay

Lead Ophthalmic Science Practitioner,
Ophthalmic Imaging Department,
Royal Hallamshire Hospital,
Sheffield (United Kingdom)






Ayan Abdullahi

Ophthalmic Science Practitioner

Madyan Qureshi


VR Clinical Fellow,
Sheffield Teaching Hospitals



Case 2:

Posterior Polar Cataract – Capsular Involvement



Patient was assessed in the cataract preoperative clinic with bilateral posterior polar cataracts. CASIA2 OCT was performed to assess the integrity of the posterior capsule. The scan revealed significant capsular involvement. This information allowed better surgical planning, and the patient had uncomplicated cataract surgery on a VR list.



Case 3:

Subluxed Lens and Vitreous Prolapse Following Trauma

The patient had blunt force trauma resulting in traumatic mydriasis, zonular trauma, subluxed lens, and vitreous prolapse into the anterior chamber. The CASIA2 OCT demonstrated the relative position of the lens, iris, and vitreous prolapse.





Dr. Sophia Reifeltshammer, MD
Department of Ophthalmology,
Kepler University Clinic, Linz (Austria)



Dr. Klaus Straßmair, MD
Department of Ophthalmology,
Kepler University Clinic, Linz (Austria)

Case 1:

Prolonged Toxic Anterior Segment Syndrome

A 74-year-old female presented with acute visual deterioration in the left eye 10 days after uncomplicated cataract surgery. Visual acuity was positive finger counting, with an intraocular pressure of 17 mmHg. Slit-lamp examination showed a severe inflammatory reaction of the anterior chamber with cells SUN 3+ and a fibrin membrane in the pupillary plane, which limited fundus view as well as the view on the implanted IOL. Ultrasound examination showed no signs of endophthalmitis. CASIA2 was used to document and visualise the extent of the inflammation. The fibrin layer was clearly visualised through the 18 radial sections and showed an attachment of the membrane to the paracentesis. The severe inflammatory reaction in the anterior chamber could also be visualised.

She was then hospitalised with suspected prolonged toxic anterior segment syndrome and received topical corticosteroids as well as pressure regulators. The progression was continuously monitored using the CASIA2 and showed improvement throughout the stay. The OCT images of the anterior chamber inflammatory reaction clearly correlated with the clinical findings and was used to objectively assess the response to medical treatment.




Figure 1: Day 1 assessment




Figure 2: Day 3 assessment




Figure 3: Day 7 assessment with a resolved membrane

Case 2:

Iris Incarceration Following Difficult Cataract Surgery

A 65-year-old female presented for a routine visit to our clinic three weeks after surgery. Cataract surgery was slightly more difficult due to a shallow anterior chamber and narrow pupil. The patient was satisfied after surgery and claimed to be free of symptoms. During the slit-lamp examination, an iris incarceration was observed at the paracenteses at the 10 o'clock position. A CASIA2 anterior segment OCT was performed for a more detailed assessment. The radial cross-sectional images enabled a clear visualisation of the incarceration and the underlying structures. The imaging allowed precise planning of the anterior chamber revision on the following day. In the postoperative follow-up anterior segment OCT, a mild persistent iris defect was seen at the affected location, which was not perceived as bothersome by the patient.




Figure 4: AS-OCT visualising the iris incarceration





Figure 5: AS-OCT visualising the postoperative defect

Case 3:



Re-rotation of a Toric IOL

At our clinic, we routinely use the CASIA2 for preoperative evaluation of corneal astigmatism and subsequent planning of both toric and multifocal IOLs. Beyond its use in preoperative planning, it also provides valuable insights for the objective assessment of lens tilt and decentration before and after surgery.

A 65-year-old female patient underwent uncomplicated cataract surgery and received a toric intraocular lens to correct her pre-existing corneal astigmatism. The IOL was aligned as planned at an implantation axis of 172°.



Keratometric		
Ks	43.17 D	@ 173°
Kf	41.96 D	@ 83°
CYL	1.21 D	@ 173°
AvgK	42.56 D	
AA (10mm)	91.5 %	
ACCP	42.51 D	
Ecc (9.0mm)	0.53	
Posterior		
Ks	-6.01 D	@ 116°
Kf	-5.93 D	@ 26°
CYL	0.08 D	@ 26°
AvgK	-5.97 D	
AA (10mm)	90.4 %	
ACCP	41.51 D	
Ecc (9.0mm)	0.88	
Real		
Ks	42.25 D	@ 174°
Kf	40.87 D	@ 84°
CYL	1.38 D	@ 174°
AvgK	41.56 D	
AA (10mm)	90.4 %	
ACCP	41.51 D	
Pachymetry		
Apex	526 um	
Thinnest	(-0.85, -0.35) mm	
Pupil D	4.27 mm	
	(-0.71, -0.27) mm	
ACD [Endo.]	1.90 mm	
ESI	0 %	



Figures 6, 7, 8: Preoperative assessment of corneal astigmatism as well as the phakic IOL

At the follow-up examination one week after surgery, the patient had a residual astigmatism of -3.50 D. The axis of the toric IOL measured at the slit lamp was 130° instead of the planned implantation axis of 172°. A re-rotation of the IOL was performed on the same day without complications. The IOL then remained stable in the correct position.

The CASIA2 anterior segment OCT was used for objective assessment of the axis of the tilt and decentration of the IOL in the capsular bag before and after re-rotation.




Figure 9: Evaluation of the IOL 1 week after surgery before re-rotation



Figure 10: Evaluation of the IOL 4 weeks after the revision




Ignacio Rodriguez-Una MD, PhD

Fernandez-Vega Ophthalmological Institute, Ophthalmological Research Foundation (Fundacion de Investigacion Oftalmologica), Univerisity of Oviedo, Oviedo (Spain)


Anterior Segment Optical Coherence Tomography has been demonstrated to be a useful tool in the postoperative assessment of multiple glaucoma surgeries. This may be particularly relevant in bleb-forming techniques and in those procedures dependent of implants or other devices. In this user experience, several cases of different glaucoma surgeries are presented showing the improved imaging performance of CASIA2 and its utility.

Case 1:


Preserflo® Microshunt Glaucoma Surgery (Post-op 4 weeks)*



Anterior segment images



CASIA2 cross-sectional image using "bleb" mode



CASIA2 anterior chamber cross-sectional image

With CASIA2 it is possible to identify the presence of fluid spreading out through the tissues, leading to multilayered conjunctiva. This sign has been reported to be a surgical success factor when it is detected in the early postoperative period.


In terms of safety, one of the key points of this device is its relation with the anterior chamber structures after implantation, in particular with the corneal endothelium. CASIA2's calipre tool can be used to measure the distance from the tube to the endothelium and to the iris, as well as the total length of the tube inside the anterior chamber (from its exit throughout the irido-corneal angle).

**Santen Pharmaceutical Co., Osaka, Japan*

Case 2A:

Non-penetrating Deep Sclerectomy (NPDS) (Post-op 3 months)

Several non-absorbable implants for NPDS have been described. Due to the high-quality imaging performance and resolution of CASIA2, it is possible to clearly identify both plates of the ESNOPER CLIP®* implant: the upper plate positioned on the scleral bed and the notched plate inserted in the suprachoroidal space (red arrow heads).



Intraoperative and anterior segment images




CASIA2 cross-sectional image using "bleb" mode

Case 2B:

Non-penetrating Deep Sclerectomy (Post-op 1 week)

Differently, the ESNOPER V2000®* is a single-plate design scleral implant that can be sutured with 10-0 nylon above the scleral bed. CASIA2 provides accurate information about the scleral lake and the drainage of aqueous humor.



Intraoperative and anterior segment images



CASIA2 cross-sectional image using "bleb" mode

*AJL ophthalmic, Vitoria, Spain

Case 3:

XEN® Gel Stent Glaucoma Surgery (Post-op 2 years)*

After a subconjunctival bleb-forming surgery, it is important to assess the exact location of the implanted device and its relation to the different ocular structures. Ensuring positioning of the distal end of the XEN® Gel Stent can predict surgical outcomes. Both eyes of the same patient are presented, with different XEN® placements and bleb morphology, but with successful results in both cases.

Right eye (OD) – Sub-Tenon's XEN® Gel Stent implantation

Although the distal end of the device was not visible by the means of the slit-lamp examination, CASIA2 was able to confirm its sub-Tenon's placement (white arrow head) and the presence of filtration at that level (red arrow heads).




Anterior segment images




CASIA2 cross-sectional image using "bleb" mode

Left eye (OS) – Subconjunctival XEN® Gel Stent implantation

Positioning the distal end of the device (red arrow head) beneath the conjunctiva led to subconjunctival filtration (yellow arrow heads), and a hyporeflective bleb with microcysts. All of these features can be evaluated with CASIA2.




Anterior segment images



CASIA2 cross-sectional image using "bleb" mode

**AbbVie Inc., Chicago, Illinois, USA*



+ Glaucoma application

With our STAR 360° application, the CASIA2 automatically measures the anterior chamber angles all around the anterior segment – thanks to its automatic scleral spur detection. It provides you with specific parameters needed to detect and treat glaucoma. With the added function “Narrow Angle Index”¹⁾ you immediately receive data about a possible narrow angle plus a referring index based on normative data.

¹⁾ available from version 60, optional feature “STAR Analysis software”

Case 1:

Post-DSAEK Complications and Cataract Surgery in Fuchs' Dystrophy

This concerns a patient with bilateral Fuchs' dystrophy whose history includes DSAEK surgery in the right eye with subsequent ocular hypertension due to misdirection syndrome, requiring pars plana Vitrectomy and Goniosynechialysis in the same procedure. The patient presented significant photophobia and cataract, with intraocular pressure of 18 mmHg (figure 1).



Dr. Pau Romera

Hospital Universitari Germans Trias i Pujol,
Barcelona (Spain)




Figure 1




Figure 2

CASIA2 images inform us about corneal transparency and good adaptation of the endothelial graft, but also allows us to visualise the integrity of the lens after inferior iridectomy and gives us a 360° visualisation of the iridocorneal angle (figure 2).




Figure 3

In this case, a narrow angle with no presence of previous synechiae in all four quadrants is observed (figure 3). Given this situation, cataract surgery will be proposed with pupillary sphincter reconstruction using iris cerclage, with no associated glaucoma surgery as part of the procedure at this time.

Case 2:

Advanced Imaging of Preserflo Microshunt in Uveitic Glaucoma Treatment

This is a 50-year-old male patient who had a history of uveitic glaucoma in his left eye, with previous Non-Penetrating Deep Sclerectomy surgery. It was decided to re-intervene in this eye after presenting pressures of 30 mmHg, performing a Preserflo Microshunt surgery (Santen, Osaka, Japan) in the superior temporal quadrant.




Figure 1




Figure 2

CASIA2 images allow us to locate the intraocular portion of the implant, calculate the distance from the tip of the tube to endothelium, as well as visualise the entire tube pathway, its internal lumen (figure 1), and the presence of posterior conjunctival bleb, visualised as a hyporeflective area in the most posterior part of the implant (asterisk in figure 2).




Figure 3

We can also make a 3D reconstruction of anterior segment localising the tube entrance (figure 3). The current pressure is 10 mmHg.

Case 3:

Long-term Success of Esnoper Clip Sclerectomy in Chronic Glaucoma

This is an 80-year-old woman with chronic open-angle glaucoma with poor pressure control despite maximum pharmacological treatment. On examination, pressures of 22 mmHg are found in both eyes with terminal glaucomatous optic nerve damage. It was decided to perform Non-Penetrating Deep Sclerectomy with Esnoper Clip implant (AJL, Alava, Spain) in both eyes. After one year, diffuse well-vascularised blebs are present in both eyes (figure 1: right eye; figure 2: left eye).



Figure 1




Figure 2

CASIA2 allows us to obtain images of the conjunctival bleb, visualising the presence of subconjunctival cystic spaces, the integrity of the intra-scleral lake, and posterior transscleral flow (figure 3, right eye, arrow), as well as visualising the presence of supraciliary filtration allowed by this implant.



Figure 3





Figure 4

In the left eye section, supraciliary flow can be visualised (figure 4, left eye, arrow). Pressures after one year were 8 mmHg in the right eye and 5 mmHg in the left eye.



+ ICL surgery

Pre-OP ICL: With two different integrated size recommendation formulas, this application enables the doctor to confirm the correct sizing for ICL surgery. A formula for calculating the recommended ICL size for vertical fixation¹⁾ has also been added.

Post-OP ICL: The Post-OP ICL movie function²⁾ shows the dynamic vault and displays it under mydriasis and miosis conditions, as well as the relevant vault range difference. The fully automated analysis output provides objective vault assessment and supports quality control.

¹⁾ available from software release version 50.7D.

²⁾ available from version 60; Reference study: "Dynamic Assessment of Light-Induced Vaulting Changes of Implantable Collamer Lens with Central Port by Swept-Source OCT: Pilot Study" Dr. Felix Gonzalez-Lopez; Madrid, Spain.



Dr. Teresa Sandmayr, FEBO

Department of Ophthalmology,
Barmherzige Brüder,
Konventhospital, Linz (Austria)




Dr. Aron Roman Czapek

Department of Ophthalmology,
Barmherzige Brüder,
Konventhospital, Linz (Austria)


Case:

PTK for Treatment of Foreign Body-Related Scarring


This case involves a 47-year-old male patient who presented with blurred vision and glare due to a central corneal scar following the removal of a metal foreign body. We used the CASIA2 to assess the depth of the lesion and to plan the PTK procedure.




Slit lamp examination shows the corneal scar prior to PTK



Anterior OCT showing a hyperdense area in the anterior corneal stroma prior to PTK



Slit lamp examination revealed almost no residual findings after PTK



Anterior OCT showing an almost clear cornea after PTK

This case demonstrates the potential effectiveness of PTK in treating corneal scarring. The patient initially had a best corrected visual acuity (BCVA) of 0.8 in the affected eye. After undergoing PTK, the BCVA improved to 1.2, and the patient was highly satisfied with the improvement in vision.



Dr. Stephan Szegedi

Department of Ophthalmology,
Clinic Hietzing, Vienna Healthcare Group,
Vienna (Austria), Karl Landsteiner Institute
for Process Optimisation and Quality
Management in Cataract Surgery,
Vienna (Austria)

Case 1:

Iatrogenic Injury of the Lens Capsule

An 80-year-old woman presented with worsening vision five days after an intravitreal anti-vascular endothelial growth factor injection for age-related macular neovascularisation. The patient's file noted early dementia and reduced cooperation during intravitreal injections. Best-corrected visual acuity (BCVA) had declined from 0.5 pre-injection to 0.05 at presentation. Slit-lamp examination revealed a clear cornea, quiet anterior chamber, nuclear cataract, and probable defect of the posterior lens capsule (Figures 1 and 2). The fundus view was limited, but no signs of endophthalmitis were recorded.



Figure 1



Figure 2

Anterior segment OCT demonstrated a posterior lens capsule injury (Figure 3), suspected to be iatrogenic from the recent intravitreal injection. This diagnosis guided the decision to perform lens extraction under general anesthesia with vitrectomy standby.




Figure 3

Intraoperatively, the capsular defect led to nucleus drop into the vitreous cavity, necessitating vitrectomy and sulcus intraocular lens implantation. Postoperatively, BCVA improved to 0.63 after nine days.

This case highlights the utility of high-resolution anterior segment OCT in definitively diagnosing posterior capsule injuries. The imaging modality provided objective preoperative information, allowing for appropriate surgical planning in managing this uncommon complication of intravitreal injections.

Case 2:

Pupillary Block Caused by Phakic Intraocular Lens

A 53-year-old man was referred for evaluation following a five-day history of intermittent pressure sensation and vision loss in the left eye. The referring ophthalmologist had initially recorded an intraocular pressure (IOP) of 35 mmHg OS and initiated treatment with oral acetazolamide. Twenty-two years prior, at age 31, the patient had undergone bilateral implantation of angle-supported phakic intraocular lenses (ACIOLs) with superior iridectomies to address high myopia.

Upon examination, best-corrected visual acuity was 0.4 OD and 0.63 OS. IOP had normalised to 10 mmHg OD and 12 mmHg OS with the acetazolamide treatment. Slit-lamp examination revealed a shallow anterior chamber (AC) in the left eye with iris bombé configuration. The iris also showed signs of pigment dispersion. The superior iridectomies performed during the original ACIOL implantation were not visible.

Anterior segment OCT demonstrated angle closure in approximately one-third of the circumference in the left eye, confirming iris bombé in the inferonasal and superotemporal quadrants (Figure 1). The right eye showed an open-angle configuration (Figure 2). The OCT also visualised the ACIOL position and its relationship to the iris. Gonioscopy confirmed these findings.




Figure 1




Figure 2

A fundus examination and optic nerve head OCT revealed borderline glaucomatous changes in both eyes, more pronounced in the left.

Based on these findings, surgical intervention was planned. The anterior chamber IOL was explanted through a scleral tunnel approach. This was followed by an uneventful phacoemulsification with implantation of a posterior chamber IOL in the capsular bag.

This case highlights the utility of anterior segment OCT in objectively documenting and quantifying the extent of angle closure in late-onset pupillary block following ACIOL implantation, even in the presence of a surgical iridectomy and normalised IOP after initial treatment. OCT imaging was the ideal means for providing documentation and objective information for treatment planning. The case also underscores the importance of long-term follow-up in patients with ACIOLs.

Case 3:

Management of Capsular Bag Phimosis

A 63-year-old woman was referred for laser capsulotomy for capsular phimosis diagnosed three months after uneventful bilateral cataract surgery. During examination, pinhole capsular phimosis was noted.

Anterior segment OCT imaging was employed to document and visualise the extent and unusual thickness of the phimotic capsular bag.

Nd:YAG laser capsulotomy was performed over three sessions, with 150-600 mJ of laser energy applied per eye per session. OCT was utilised to monitor the progress of the interrupted surgical treatment, allowing for documentation of the capsulotomy's extent and effectiveness (Figure 1: after first capsulotomy, Figure 2: follow-up visit after final treatment).



Figure 1




Figure 2

For postoperative treatment, NSAID and cortisone eye drops were prescribed. Despite this, the patient developed mild cystoid macular edema postoperatively. With continued treatment, macular edema resolved over a one-month period.

OCT was utilised to document the extent and unusual thickness of the phimotic capsular bag and allowed for assessment of the progress during the multi-session treatment.

The background of the page is a soft-focus photograph of a mountain range. The mountains are layered, with the closest peaks in shades of blue and the more distant ones fading into a pale, hazy sky. The overall tone is serene and atmospheric.

Join this project

*Do you also have interesting cases that could
be valuable for other ophthalmologists?
Then feel free to contribute to this brochure.*

*Please contact us at **marketing@tomey.de** –
we really look forward to hearing about your
practical experiences with our AS-OCT CASIA2.*




Ana Trpeska Boshkoska MD
Ophthalmologist, PHI University Clinic
for Eye Diseases, Skopje (N. Macedonia),
Teaching assistant

Case 1:

Peters Anomaly

This CASIA2 OCT tomographic image shows a case of a 58-year-old patient with unilateral Peters Anomaly. Corneal opacity and iridocorneal adhesions are present. There are no lenticulocorneal adhesions.




Topography

ID:

Name:

TOMEY
CASIA2 Ver.50.5A.03

Exam Date:
Comment:




Clinic Info.: Univerzitetna Klinika za Očni bolesti / Skopje

Case 2:

Phacolytic Glaucoma


With CASIA2, the presence of lens proteins and protein-containing macrophages in the anterior chamber is easily noticeable, as in this case of a patient with a phacolytic glaucoma. Hypermature cataract is visible. Corneal edema is not yet established. There is also advanced pterygium verum present.



Case 3:

AC-IOL and Peripheral Iridectomy Not Functioning

This is a case of a patient that underwent cataract surgery and AS-IOL was implanted. The peripheral iridectomy is not functioning, iris bombe is present, intraocular pressure is raised and the AS-IOL is pushed anteriorly. Corneal edema is established.





Dr. Jorge Vila Arteaga
La Fe Hospital, Valencia (Spain)

Case 1:

Descemet's Detachment – Rare Post-NPGS Complication

A 65-year-old patient underwent deep sclerectomy – non-penetrating glaucoma surgery (NPGS) – for glaucoma five years ago. The IOP was maintained at levels between 11-13 mmHg during this time. He went to the emergency room due to decreased visual acuity in the operated eye. IOP was 12 mmHg.



Figure 1

Corneal edema was observed and treatment was started with topical corticosteroids and hyperosmolar agents. He was referred to the Glaucoma Unit.

On examination we found a superior, semicircular, and well-differentiated corneal edema (Figure 1). For the correct diagnosis, AS-OCT is performed using CASIA2.

In the surgical area, the scleral lake (SL) and the trabeculo-Descemet window (TDW) can be seen. There is corneal thickening and detachment of Descemet's membrane (***) (Figure 2).




Figure 2

Descemet's membrane detachment is an unusual complication of NPGS that is often associated with an increase in IOP due to a tenon cyst. On this occasion the IOP was within normal limits. Its diagnosis is very complex if AS-OCT is not available, because with the slit lamp it is not possible to differentiate the Descemet's membrane due to the clouding of the cornea. The treatment consists of injection of 20% SF6 in the anterior chamber (Figure 3) which manages to apply the membrane. This is constable using the AS-OCT (Figure 4).



Figure 3







Figure 4

Case 2:

Post-Surgery Corneal Edema and Athalamia in Glaucoma

A 45-year-old woman underwent phakic anterior chamber lens surgery 23 years ago. She develops glaucoma and undergoes a trabeculectomy. The day after surgery the chamber is wide, the IOP is 5 mmHg.

At the weekly check-up, the patient presented with corneal edema and athalamia. She underwent intensive treatment with mydriatic agents and corticosteroids. Performing an AS-OCT using CASIA2 allows all parameters to be monitored qualitatively and quantitatively: corneal edema, distance from the lens to the endothelium(*), irido-corneal contact, and the inflammatory pupillary membrane(**).



After 30 minutes the improvement was already visible on the OCT, and the next day the improvement could be quantified. In 3 days the clinical condition was resolved.

Case 1:

Macular Corneal Dystrophy



A 45-year-old female patient with a significant family history of macular dystrophy presented with complaints of reduced vision. Her best corrected visual acuity (BCVA) was 0.8 in the right eye (OD) and 0.2 in the left eye (OS).

The right eye had undergone penetrating keratoplasty (PKP) 7 months prior due to corneal pathology. The left eye had no surgical history but presented with progressive visual deterioration.






Prof. Adam Wylęgała

Experimental Ophthalmology Unit,
Department of Biophysics,
Medical University of Silesia,
Railway Hospital, Katowice (Poland)



The left eye showed generalised corneal thinning and the increased corneal stromal reflectivity characteristic of macular dystrophy.






The graft appeared structurally intact with acceptable thickness centrally and regular posterior corneal contour. AvgK post-PKP was 41.49 D, with residual astigmatism of 2.79 D. The CASIA2 anterior segment OCT proved crucial in evaluating corneal structure and detecting macular corneal dystrophy, particularly in the left eye. The right eye showed a stable post-PKP graft with residual regular astigmatism.

Case 2:




Pellucid Marginal Degeneration

A 74-year-old female patient presented with gradually decreasing vision, particularly in her right eye, with notable difficulties at night and an increase in astigmatism that her regular glasses no longer corrected adequately. The best corrected visual acuity (BCVA) in the right eye (OD) was 0.2 with high against-the-rule astigmatism, approximately $-12.00D \times 110^\circ$, and in the left eye (OS) it was 0.5 with moderate against-the-rule astigmatism.



Topographic examination revealed a classic “crab claw” pattern on axial power maps, particularly in the right eye, with marked inferior corneal steepening and flattening in the vertical meridian above the thinned area. Against-the-rule astigmatism was evident, with the right eye showing a steep K (Ks) of approximately 51D and a flat K (Kf) of around 35 D, resulting in significant corneal astigmatism (>15 D). In the left eye, the astigmatism was less pronounced but still significant (6-8 D). The minimum corneal thickness was approximately 450 μm in the right eye and 500 μm in the left eye.

Vertical scans displayed normal central corneal architecture with abrupt thinning in the inferior peripheral cornea. CASIA2 imaging confirmed the diagnosis of pellucid marginal degeneration rather than inferior keratoconus, characterised by peripheral rather than paracentral thinning and a distinctive crab claw topographic pattern.



This trend analysis compares keratometric data over a span of 4 years and 1 month.

The average keratometric power has remained consistent, with minimal variation of 0.87 D in the right and 0.08 D in the left eye.



Azamat Farkhadovich Yusupov, MD, DSc, Professor

Director of the Republican Specialised Scientific and Practical Medical Centre of Eye Microsurgery, Chief Ophthalmologist of the Ministry of Health of the Republic of Uzbekistan



Ikbol Fazilovich Saliev, MD, PhD

Author and originator of the surgical concept of dual intraocular lens implantation, Founder of SAIF Optima Eye Clinic, Tashkent (Uzbekistan), distinguished ophthalmic surgeon and mentor (in memoriam)



Farzona Abdurakhimova, MD

Ophthalmic surgeon and researcher, SAIF Optima Eye Clinic, Tashkent (Uzbekistan)

Case:

Visualisation of Dual Intraocular Lenses Using CASIA2 Anterior Segment OCT

A patient with high myopia underwent refractive lens exchange with implantation of two intraocular lenses within the capsular bag. Dual intraocular lens implantation may be used in selected clinical cases to achieve the required refractive correction in patients with high myopia and extremely long axial length.

However, during routine slit-lamp examination, it is often difficult to clearly visualise the spatial configuration of two intraocular lenses located inside the capsular bag. To evaluate the anatomical structures of the anterior segment more precisely, imaging was performed using CASIA2 anterior segment OCT.

CASIA2 Examination

Anterior segment OCT imaging provided high-resolution cross-sectional images of the ocular structures.

The use of CASIA2 allowed clear visualisation of:

- + the cornea
- + the anterior chamber
- + the capsular bag
- + the position of the two intraocular lenses

The OCT images clearly demonstrated the configuration of the two intraocular lenses positioned within the capsular bag.

Measurements

CASIA2 also allowed accurate measurements of important anatomical parameters of the anterior segment, including:

- + anterior chamber depth
- + distance from the cornea to the posterior capsule
- + spatial configuration of the intraocular lens complex

These measurements provide objective information about the position of the lenses after surgery.

Discussion

Dual intraocular lens implantation can be an effective method for refractive correction in highly myopic eyes. However, evaluation of the spatial configuration of two intraocular lenses inside the capsular bag can be challenging using conventional slit-lamp examination. Anterior segment OCT imaging provides high-resolution cross-sectional visualisation of the anterior segment and allows precise assessment of intraocular lens position and configuration.

Accurate visualisation of the intraocular lens complex is important for post-operative evaluation and for early detection of possible complications such as IOL decentration or interlenticular changes.




Figure 1: Intraoperative microscopic image demonstrating the configuration of two intraocular lenses positioned within the capsular bag during dual IOL implantation.



Figure 2: Anterior segment OCT image obtained with CASIA2 demonstrating the position of dual intraocular lenses within the capsular bag and measurement of the distance from the cornea to the posterior capsule.



Figure 3: CASIA2 anterior segment OCT lens analysis demonstrating the configuration, centration, and spatial position of dual intraocular lenses within the capsular bag.



Figure 4: Comparison of anterior segment OCT images obtained with CASIA2 demonstrating the anatomical configuration of the anterior segment and the position of the intraocular lens complex.

Clinical Value

Conventional diagnostic methods do not always allow precise visualisation of the spatial relationship between intraocular structures. The use of CASIA2 anterior segment OCT provides detailed imaging of the anterior segment and allows clear demonstration of the configuration of two intraocular lenses within the capsular bag. With this technology it is possible to measure the distance from the cornea to the posterior capsule and objectively evaluate the anatomical position of the intraocular lenses.

Educational Value

High-quality CASIA2 images are also extremely useful for scientific presentations and educational purposes. During scientific meetings it is not always possible to demonstrate surgical videos due to time limitations of presentations. Anterior segment OCT images allow surgeons to clearly demonstrate the position of intraocular lenses and explain the surgical concept quickly and effectively.

Conclusion

The use of CASIA2 anterior segment OCT allows precise visualisation of dual intraocular lenses within the capsular bag and enables accurate measurements of anterior segment anatomy. This technology represents a valuable diagnostic and educational tool for documenting the results of complex intraocular surgical procedures.

Acknowledgement

The author expresses sincere gratitude to her late mentor Ikbol Fazilovich Saliev, who first developed and patented the surgical concept of dual intraocular lens implantation. His innovative ideas laid the foundation for the surgical approach presented in this report.


+ Artificial prosthetic eyes

apEye


The CASIA2 apEye software enables a fully digital workflow for the 3D production of ocular prostheses. The process is simple, efficient and pain-free for patients. Doctors, ophthalmologists, hospitals and clinics also benefit from a streamlined workflow and the latest advances in digital manufacturing. The result is a 3D printed ocular prosthesis with true-to-life colours and a unique iris topography.

Partners

The apEye project brought together leading partners from research, medicine and additive manufacturing. Software development was supported by Fraunhofer IGD in Darmstadt, while medical expertise was provided by Moorfields Eye Hospital in London. Together, they covered the entire workflow – from image acquisition and data processing to 3D printing and clinical evaluation. Today, the project is continued by Kröner Medical under the name ocla. Kröner Medical is responsible for the distribution of the ocular prostheses and the coordination of production, logistics and shipping.



Fellow eye scan (AS Global scan)



TOMEY EUROPE
TOMEY GMBH

Wiesbadener Strasse 21
90427 Nuremberg | Germany
+49 911 938 546 2 - 0
info@tomey.de

tomey.de



Follow TOMEY



The statements of the ophthalmologists who provided testimonials only reflect their personal opinions and experiences and do not necessarily reflect the opinions with whom they are affiliated. The ophthalmologists giving these testimonials may have a contractual relationship with TOMEY GmbH and may have received financial compensation.

Always read and follow the instructions for use.
Not all products, services or offers are approved or offered in every market.
Please note that the current status of approval for the labelling, instructions and contents of the brochure may vary from one country to another.

TOMEY GmbH is the European headquarters of TOMEY Corporation,
2-11-33 Noritakeshinmachi Nishi-Ku, Nagoya, 451-0051, Japan